



TOXICOLOGY POCT DEVICES

FROM YOUR LABORATORY SERVICES PROVIDER

PAML offers Point of Collection Testing (POCT) devices to detect drugs of abuse in human urine. For Workplace Drug Testing an Integrated Split Key Cup is available with two different configurations.

Products available are:

Name: Split Specimen 6 Panel Cup w/Adulteration
Drugs: (AMP/THC/COC/mAMP/MOR/PCP)
PAML Part #: 8355

Name: Split Specimen 10 Panel Cup w/Adulteration
Drugs: (mAMP/COC/THC/AMP/OPI/OXY/PCP/BZO/BAR/MTD)
PAML Part #: 8354



Drug Name:	Abbreviation:	Cutoff Levels:
Amphetamine	AMP	1000 ng/mL
Barbiturates	BAR	300 ng/mL
Benzodiazepines	BZO	300 ng/mL
Cocaine	COC	300 ng/mL
Marijuana	THC	50 ng/mL
Methadone	MTD	300 ng/mL
Methamphetamine	mAMP	1000 ng/mL
Morphine	MOR	300 ng/mL
Opiates	OPI	2000 ng/mL
Oxycodone	OXY	100 ng/mL
Phencyclidine	PCP	25 ng/mL

Provided for the clients of

Providence Associates Medical Laboratories PacLab Network Laboratories
Tri-Cities Laboratory Treasure Valley Laboratory Alpha Medical Laboratory
For more information, please contact your local representative.

Features:

- Integrated collection cup and test cassette allows for easy urine collection and prevents tampering.
- Each cup has a temperature strip affixed to its side to enable users to verify the validity of the specimen as the urine is provided.
- Test results are kept confidential until the authorized individual is ready to view them.
- Collectors may record pertinent information regarding the specimen directly on the collection cup to aid in identification.
- Custom designed test cup provides for easy handling during the collection process.
- Split specimen for the testing procedure ensures remaining specimen is unaffected in case it is needed for confirmatory testing.
- Test devices are self-contained and portable and can be sent directly to the laboratory if confirmation is required.
- Test controls are built into the device.
- Negative results are available in as few as 3 minutes, with positive results in 5 minutes.
- Results appear clear and are easy-to-read.

Easy To Use

1 Collect Specimen



2 Insert Key




3 Peel Back Label



4 Read the Results



Read results at (5) minutes. Results stable up to (60) minutes.



Adulteration Strip Identifies Oxidants, Specific Gravity, and pH



For more information, please call PAML Toxicology Client Services at 1-800-541-7891 x7 or (509) 755-8991.